

Magog, 9 May 2018



SliceOmatic balances powerful segmentation tools with a highly interactive interface, enabling you to completely control the segmentation process. Different modalities or even different tissues within the same slice are best segmented with different techniques. SliceOmatic makes it easy to use the appropriate technique for each specific task.

The process of identifying each tissue inside an image is known as "segmentation". SliceOmatic has highly accurate segmentations of the muscle, SAT (Subcutaneous Adipose Tissue) and VAT (Visceral Adipose Tissue) regions.

SliceOmatic is only available from the manufacturer, Tomovision.



Yves Martel
President