

# PULSE

PU-C-AA (shown)



modular

## Pulse *(see online price list for complete specs)*

CODE	ITEM	W	D	AH	BH	SW	SD	SH	WT
PU-C-RA	Chair	30.5	31	24	34	21	20.5	19	84
PU-L-RA	Loveseat	51.5	31	24	34	42	20.5	19	94
PU-S-RA	Sofa	72.5	31	24	34	63	20.5	19	111

Bariatric Chairs and Modular Configurations Available in all Arm Styles

W=Width D=Depth AH=Arm Height BH=Back Height SW=Seat Width SD=Seat Depth SH=Seat Height WT=Weight

## PULSE SPECIFICATIONS **signature series**

### BACK

The back frame is constructed with 7/8" engineered paneling with the support system constructed of 11-gauge sinuous wire linked together with a 12-gauge tie wire. The back support system is covered with woven versare fabric with one piece of high-density polymeric foam. The outside back of the frame is covered with woven tyar and an additional piece of high-density polymeric foam.

### BASE

The base frame is constructed completely of 7/8" engineered paneling. Legs are attached with steel bolts and feet have plastic or metal glides.

### SEAT

Seat cushions are of high-density, polymeric foam. The top and bottom is covered with a laminated, polyurethane, moisture-barrier membrane. The cushion is secured to the base frame by a nylon security strap (steel cable optional). Uniform seating comfort is provided by a dual-leverage amplification Spring System.

### ARM

Arm frames are constructed completely of 7/8" engineered paneling and 80-point chipboard. The arm frames are covered with polymeric high-density urethane foam. The arm frame is attached to the base frame by an interlocking design.

### UPHOLSTERY

All upholstery offerings conform to California Technical Bulletin 117 fire code standard. All estimated yardage quotes are based on 54" non-matched, solid color material unless otherwise specified. Additional yardage and matching fees may apply.

### ASSEMBLY AND DISASSEMBLY

No tools are needed. All components can be disassembled and removed by loosening the thumbscrew and sliding the components out of the tamper-proof interlocking state. Assembly can be completed within minutes without the usage of tools, and all upholstery fastens with hook and loop tape, facilitating removal and replacement of covers.

### WARRANTY

This product is warranted against defects in materials and workmanship and against structural damage occurred through normal use for 35-years from date of purchase. Foam is also covered under the 35-year period while fabric is generally covered under the fabric company's warranty. Transformations does not warrant any materials beyond its control.

### ADDITIONAL OPTIONS

CAL 133 fabric treatment may also be applied for an additional fee per yard.


*This is an abbreviated specifications list. For complete specs, contact Customer Service.*

### AVAILABLE OPTIONS:

Chairs (5 Arm Styles)  
Loveseats in all styles  
Sofas in all styles  
Modular Configurations  
Bariatric Chairs

### UPGRADE OPTIONS:

CAL133  
Tablet Arm (Available  
on FC and SA Arm  
Styles Only)



**SEE PRICE LIST**

**FOR MORE INFORMATION**

### FOOT OPTIONS:

STAIN CHOICES:  
Black, Hickory,  
Natural, Sable,  
Sienna, Umber

TF08



Visit our  
[Options](#)  
web page for  
samples.



**PU-C-RA**  
chair - round arm



**BARIATRIC CHAIRS**  
shown: PU-BC-SA  
available in all arm styles



**PU-L-AA**  
loveseat available in all styles

**PU-C-SA**  
chair - straight arm



**PU-C-AA**  
chair - arched arm



**PU-S-AA**  
sofa - available in all styles

**PU-C-FC**  
chair - flat cap arm



For More Information:

Transformations Furniture 310 Racquet Drive Fort Wayne IN 46825

www.trfurniture.com customerservice@trfurniture.com phone: **800.440.9337**